

## 小鼠抗 PRMT5 单克隆抗体

中文名称： 小鼠抗 PRMT5 单克隆抗体

英文名称： Anti-PRMT5 mouse monoclonal antibody

别名： protein arginine methyltransferase 5; JBP1; SKB1; IBP72; SKB1Hs; HRMT1L5

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： PRMT5

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes an enzyme that belongs to the methyltransferase family. The encoded protein catalyzes the transfer of methyl groups to the amino acid arginine, in target proteins that include histones, transcriptional elongation factors and the tumor suppressor p53. This gene plays a role in several cellular processes, including transcriptional regulation, and the assembly of small nuclear ribonucleoproteins. A pseudogene of this gene has been defined on chromosome 4. Alternative splicing results in multiple transcript variants encoding different isoforms.
<b>Applications:</b>	WB

<b>Name of antibody:</b>	PRMT5
<b>Immunogen:</b>	Fusion protein of human PRMT5
<b>Full name:</b>	protein arginine methyltransferase 5
<b>Synonyms:</b>	JBP1; SKB1; IBP72; SKB1Hs; HRMT1L5
<b>SwissProt:</b>	O14744
<b>WB Predicted band size:</b>	73 KD
<b>WB Positive control:</b>	K562 cell lysate
<b>WB Recommended dilution:</b>	200-1000