

兔抗 TLN1 多克隆抗体

中文名称：兔抗 TLN1 多克隆抗体

英文名称：Anti-TLN1 rabbit polyclonal antibody

别名：TLN; ILWEQ

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TLN1

克隆类型：rabbit polyclonal

反应种属：Human, Mouse

标记物：Unconjugate

技术规格

Background:	This gene encodes a cytoskeletal protein that is concentrated in areas of cell-substratum and cell-cell contacts. The encoded protein plays a significant role in the assembly of actin filaments and in spreading and migration of various cell types, including fibroblasts and osteoclasts. It codistributes with integrins in the cell surface membrane in order to assist in the attachment of adherent cells to extracellular matrices and of lymphocytes to other cells.
Applications:	ELISA, WB, IHC
Name of antibody:	TLN1

Immunogen:	Synthetic peptide of human TLN1
Full name:	talin 1
Synonyms :	TLN; ILWEQ
SwissProt:	Q9Y490
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human breast cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	270 kDa
WB Positive control:	Huvec cells
WB Recommended dilution:	500-2000



