

## 小鼠抗 PRPS1L1 单克隆抗体

中文名称： 小鼠抗 PRPS1L1 单克隆抗体

英文名称： Anti-PRPS1L1 mouse monoclonal antibody

别名： phosphoribosyl pyrophosphate synthetase 1-like 1; PRPS1; PRPS3; PRPSL; PRS-III

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： PRPS1L1

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	This intronless gene is specifically expressed in the testis, and encodes a protein that is highly homologous to the two subunits of phosphoribosylpyrophosphate synthetase encoded by human X-linked genes, PRPS1 and PRPS2. These enzymes convert pyrimidine, purine or pyridine bases to the corresponding ribonucleoside monophosphates. In vitro transcription/translation and site-directed mutagenesis studies indicate that translation of this mRNA initiates exclusively at a non-AUG (ACG) codon.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	PRPS1L1
<b>Immunogen:</b>	Fusion protein of human PRPS1L1

<b>Full name:</b>	phosphoribosyl pyrophosphate synthetase 1-like 1
<b>Synonyms:</b>	PRPS1; PRPS3; PRPSL; PRS-III
<b>SwissProt:</b>	P21108
<b>WB Predicted band size:</b>	35 KD
<b>WB Positive control:</b>	293T, Raji, MCF7 cell lysates
<b>WB Recommended dilution:</b>	500-2000