

兔抗 TLR6 多克隆抗体

中文名称：兔抗 TLR6 多克隆抗体

英文名称： Anti-TLR6 rabbit polyclonal antibody

别名： CD286

相关类别： 一抗

抗原： TLR6

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a member of the Toll-like receptor (TLR) family which plays a fundamental role in pathogen recognition and activation of innate immunity. TLRs are highly conserved from Drosophila to humans and share structural and functional similarities. They recognize pathogen-associated molecular patterns (PAMPs) that are expressed on infectious agents, and mediate the production of cytokines necessary for the development of effective immunity. The various TLRs exhibit different patterns of expression. This receptor functionally interacts with toll-like receptor 2 to mediate cellular res

	ponse to bacterial lipoproteins.
Applications:	ELISA, IHC
Name of antibody:	TLR6
Immunogen:	Fusion protein of human TLR6
Full name:	toll-like receptor 6
Synonyms :	CD286
SwissProt:	Q9Y2C9
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human breast cancer
IHC Recommend dilution:	50-200

