

小鼠抗 LRRFIP1 单克隆抗体

中文名称： 小鼠抗 LRRFIP1 单克隆抗体

英文名称： Anti-LRRFIP1 mouse monoclonal antibody

别名： LRR binding FLII interacting protein 1; GCF2; TRIP; FLAP1; GCF-2; FLAP-1; HUFI-1; FLIAP1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： LRRFIP1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:	Transcriptional repressor which preferentially binds to the GC-rich consensus sequence (5'-AGCCCCGGCG-3') and may regulate expression of TNF, EGFR and PDGFA. May control smooth muscle cells proliferation following artery injury through PDGFA repression. May also bind double-stranded RNA. Positively regulates Toll-like receptor (TLR) signaling in response to agonist probably by competing with the negative FLII regulator for MYD88-binding.
Applications:	WB

Name of antibody:	LRRFIP1
Immunogen:	Fusion protein of human LRRFIP1
Full name:	LRR binding FLII interacting protein 1
Synonyms:	GCF2; TRIP; FLAP1; GCF-2; FLAP-1; HUFI-1; FLIAP1
SwissProt:	Q32MZ4
WB Predicted band size:	89 kDa
WB Positive control:	Raji, C2C12, MCF7 and Jurkat cell lysates
WB Recommended dilution:	1000-5000

