

小鼠抗 LSM11 单克隆抗体

中文名称： 小鼠抗 LSM11 单克隆抗体

英文名称： Anti-LSM11 mouse monoclonal antibody

别名： LSM11, U7 small nuclear RNA associated

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： LSM11

宿主： Mouse

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	Component of the U7 snRNP complex that is involved in the histone 3'-end pre-mRNA processing (By similarity). Increases U7 snRNA levels but not histone 3'-end pre-mRNA processing activity, when overexpressed. Required for cell cycle progression from G1 to S phases. Binds specifically to the Sm-binding site of U7 snRNA.
Applications:	WB
Name of antibody:	LSM11
Immunogen:	Fusion protein of human LSM11
Full name:	LSM11, U7 small nuclear RNA associated
SwissProt:	P83369



WB Predicted band size:	40 KD
WB Positive control:	Hela, Jurkat, HepG2 cell lysates
WB Recommended dilution:	200-1000