

兔抗 TLX1 多克隆抗体

中文名称: 兔抗 TLX1 多克隆抗体

英文名称: Anti-TLX1 rabbit polyclonal antibody

别名: T-cell leukemia homeobox 1; TCL3; HOX11

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: TLX1

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a nuclear transcription factor that belongs to the NK-linked or NK-like (NKL) subfamily of homeobox genes. The encoded protein is required for normal development of the spleen during embryogenesis. This protein is also involved in specification of neuronal cell fates. Ectopic expression of this gene due to chromosomal translocations is associated with certain T-cell acute lymphoblastic leukemias. Alternate splicing results in multiple transcript variants.
Applications:	ELISA, IHC

Name of antibody:	TLX1
Immunogen:	Synthetic peptide of human TLX1
Full name:	T-cell leukemia homeobox 1
Synonyms:	TCL3; HOX11
SwissProt:	P31314
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	10-50

