

## 小鼠抗 LZTS1 单克隆抗体

中文名称： 小鼠抗 LZTS1 单克隆抗体

英文名称： Anti-LZTS1 mouse monoclonal antibody

别名： leucine zipper tumor suppressor 1; F37; FEZ1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： LZTS1

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes a tumor suppressor protein that is ubiquitously expressed in normal tissues. In uveal melanomas, expression of this protein is silenced in rapidly metastasizing and metastatic tumor cells but has normal expression in slowly metastasizing or nonmetastasizing tumor cells. This protein may have a role in cell-cycle control by interacting with the Cdk1/cyclinB1 complex. This gene is located on chromosomal region 8p22. Loss of heterozygosity (LOH) in the 8p arm is a common characteristic of many types of cancer.
<b>Applications:</b>	WB

<b>Name of antibody:</b>	LZTS1
<b>Immunogen:</b>	Fusion protein of human LZTS1
<b>Full name:</b>	leucine zipper tumor suppressor 1
<b>Synonyms:</b>	F37; FEZ1
<b>SwissProt:</b>	Q9Y250
<b>WB Predicted band size:</b>	67 KD
<b>WB Positive control:</b>	Jurkat cell lysate
<b>WB Recommended dilution:</b>	200-1000