

## 小鼠抗 MAGEB1 单克隆抗体

中文名称: 小鼠抗 MAGEB1 单克隆抗体

英文名称: Anti-MAGEB1 mouse monoclonal antibody

别名: CT3.1; DAM10; MAGE-Xp; MAGEL1

抗原: MAGEB1

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human , Dog , Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

#### Background:

This gene is a member of the MAGEB gene family. The members of this family have their entire coding sequences located in the last exon, and the encoded proteins show 50 to 68% sequence identity to each other. The promoters and first exons of the MAGEB genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. This gene is localized in the DSS (dosage-sensitive sex reversal) critical region, and expressed in testis and in a significant fraction of tumors of various histological types. Thi

	s gene and other MAGEB members are clustered on chromosome Xp22-p21. Multiple alternatively spliced transcript variants encoding the same protein have been found for this gene, however, the full length nature of some variants has not been defined. [provided by RefSeq, Jul 2008].
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	MAGEB1
<b>Immunogen:</b>	Fusion protein of human MAGEB1
<b>Full name:</b>	melanoma antigen family B1 (MAGEB1), transcript variant 2
<b>Synonyms:</b>	CT3.1; DAM10; MAGE-Xp; MAGEL1
<b>SwissProt:</b>	P43366
<b>IHC positive control:</b>	carcinoma of human bladder tissue
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	39 kDa
<b>WB Positive control:</b>	HepG2, HeLa, SVT2, A549, Jurkat, MDCK cell lysates
<b>WB Recommended dilution:</b>	1000-5000