

兔抗 TMEM8B 多克隆抗体

中文名称：兔抗 TMEM8B 多克隆抗体

英文名称：Anti-TMEM8B rabbit polyclonal antibody

别名：NGX6, NAG-5, C9orf127, RP11-112J3.10

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TMEM8B

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May function as a regulator of the EGFR pathway. Probable tumor suppressor which may function in cell growth, proliferation and adhesion.
Applications:	ELISA, IHC
Name of antibody:	TMEM8B
Immunogen:	Synthetic peptide of human TMEM8B
Full name:	Transmembrane protein 8B
Synonyms :	NGX6, NAG-5, C9orf127, RP11-112J3.10
SwissProt:	A6NDV4
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human colon cancer and Human gastric cancer

IHC Recommend dilution:	25-100
-------------------------	--------

