

小鼠抗 MAP2K6 单克隆抗体

中文名称： 小鼠抗 MAP2K6 单克隆抗体

英文名称： Anti-MAP2K6 mouse monoclonal antibody

别名： mitogen-activated protein kinase kinase 6; MEK6; MKK6; MAPKK6; PRKMK6; SAPKK3; SAPKK-3

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： MAP2K6

反应种属： Human, Rat

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:

This gene encodes a member of the dual specificity protein kinase family, which functions as a mitogen-activated protein (MAP) kinase kinase. MAP kinases, also known as extracellular signal-regulated kinases (ERKs), act as an integration point for multiple biochemical signals. This protein phosphorylates and activates p38 MAP kinase in response to inflammatory cytokines or environmental stress. As an essential component of p38 MAP kinase mediated signal transduction pathway, this gene is involv

	ed in many cellular processes such as stress induced cell cycle arrest, transcription activation and apoptosis.
Applications:	WB, IP
Name of antibody:	MAP2K6
Immunogen:	Fusion protein of human MAP2K6
Full name:	mitogen-activated protein kinase kinase 6
Synonyms:	MEK6; MKK6; MAPKK6; PRKMK6; SAPKK3; SAPKK-3
SwissProt:	P52564
WB Predicted band size:	37 kDa
WB Positive control:	PC12, C6, Ramos, Jurkat and Hela cell lysates
WB Recommended dilution:	500-2000



