

## 兔抗 TMEM237 多克隆抗体

中文名称：兔抗 TMEM237 多克隆抗体

英文名称： Anti-TMEM237 rabbit polyclonal antibody

别名： JBTS14; ALS2CR4

抗原： TMEM237

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a tetraspanin protein that is thought to be involved in WNT signaling. Defects in this gene are a cause of Joubert syndrome-14. Two transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	TMEM237
<b>Immunogen:</b>	Synthesized peptide derived from internal of human ALS2CR4.
<b>Full name:</b>	transmembrane protein 237
<b>Synonyms :</b>	JBTS14; ALS2CR4

<b>SwissProt:</b>	Q96Q45
<b>WB Predicted band size:</b>	46 kDa
<b>WB Positive control:</b>	Jurkat cells lysate
<b>WB Recommended dilution:</b>	500-3000

