

兔抗 TMEM254 多克隆抗体

中文名称：兔抗 TMEM254 多克隆抗体

英文名称：Anti-TMEM254 rabbit polyclonal antibody

别名：C10orf57; bA369J21.6

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TMEM254

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein belongs to the TMEM214 family. Subcellular location is Membrane; Multi-pass membrane protein.
Applications:	ELISA, WB
Name of antibody:	TMEM254
Immunogen:	Synthetic peptide of human TMEM254
Full name:	transmembrane protein 254
Synonyms:	C10orf57; bA369J21.6
SwissProt:	Q8TBM7
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	14 kDa

WB Positive control:	Human heart tissue ,PC-3 cell ,Mouse heart tissue and Human cerebrum tissue lysates
WB Recommended dilution:	200-1000

