

## 小鼠抗 MED15 单克隆抗体

中文名称： 小鼠抗 MED15 单克隆抗体

英文名称： Anti-MED15 mouse monoclonal antibody

别名： mediator complex subunit 15; TIG1; CAG7A; CTG7A; PCQAP; TIG-1; TNRC7; ARC105

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： MED15

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a subunit of the multiprotein complexes PC2 and ARC/DRIP and may function as a transcriptional coactivator in RNA polymerase II transcription. This gene contains stretches of trinucleotide repeats and is located in the chromosome 22 region which is deleted in DiGeorge syndrome. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	MED15

<b>Immunogen:</b>	Fusion protein of human MED15
<b>Full name:</b>	mediator complex subunit 15
<b>Synonyms:</b>	TIG1; CAG7A; CTG7A; PCQAP; TIG-1; TNRC7; ARC105
<b>SwissProt:</b>	Q96RN5
<b>WB Predicted band size:</b>	87 KD
<b>WB Positive control:</b>	MCF7, K562 and Hela cell lysates
<b>WB Recommended dilution:</b>	500-2000