

小鼠抗 RAB3IL1 单克隆抗体

中文名称: 小鼠抗 RAB3IL1 单克隆抗体

英文名称: Anti-RAB3IL1 mouse monoclonal antibody

别名: GRAB

储存: 冷冻 (-20℃) 避光

宿主: Mouse

抗原: RAB3IL1

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	This gene encodes a guanine nucleotide exchange factor for the ras-related protein Rab3A. The encoded protein binds Rab3a and the inositol hexakisphosphate kinase InsP6K1. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 7.
Applications:	WB, FC
Name of antibody:	RAB3IL1
Immunogen:	Fusion protein of human RAB3IL1
Full name:	RAB3A interacting protein (rabin3)-like 1 (RAB3IL1), transcript variant 1
Synonyms:	GRAB

SwissProt:	Q8TBN0
WB Predicted band size:	43 kDa
WB Positive control:	SVT2, MDCK, PC12 cell lysates
WB Recommended dilution:	200-1000