

## 小鼠抗 RAB4A 单克隆抗体

中文名称： 小鼠抗 RAB4A 单克隆抗体

英文名称： Anti-RAB4A mouse monoclonal antibody

别名： RAB4A, member RAS oncogene family; RAB4; HRES1; HRES-1; HRES-1/RAB4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： RAB4A

反应种属： Human, Mouse, Rat, Monkey

标记物： Unconjugate

克隆类型： Mouse Monoclonal

### 技术规格

<b>Background:</b>	This gene is a member of the largest group in the Ras superfamily of small GTPases, which regulate membrane trafficking. The encoded protein is associated with early endosomes and is involved in their sorting and recycling. The protein also plays a role in regulating the recycling of receptors from endosomes to the plasma membrane. Alternatively spliced transcript variants have been observed for this gene.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	RAB4A
<b>Immunogen:</b>	Fusion protein of human RAB4A

<b>Full name:</b>	RAB4A, member RAS oncogene family
<b>Synonyms:</b>	RAB4; HRES1; HRES-1; HRES-1/RAB4
<b>SwissProt:</b>	P20338
<b>WB Predicted band size:</b>	24 kDa
<b>WB Positive control:</b>	COS7 cell, mouse brain tissue, Hela and PC-12 cell lysates
<b>WB Recommended dilution:</b>	500-2000

