

兔抗 TMPRSS5 多克隆抗体

中文名称：兔抗 TMPRSS5 多克隆抗体

英文名称：Anti-TMPRSS5 rabbit polyclonal antibody

别名：transmembrane serine protease 5; SPINESIN

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TMPRSS5

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein that belongs to the serine protease family. Serine proteases are known to be involved in many physiological and pathological processes. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	TMPRSS5
Immunogen:	Fusion protein of human TMPRSS5
Full name:	transmembrane serine protease 5
Synonyms:	SPINESIN
SwissProt:	Q9H3S3
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human prostate cancer and Human breast cancer
IHC Recommend dilution:	200-400



