

小鼠抗 RAB21 单克隆抗体

中文名称: 小鼠抗 RAB21 单克隆抗体

英文名称: Anti-RAB21 mouse monoclonal antibody

别名: KIAA0118

抗原: RAB21

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human , Dog , Rat , Monkey , Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	This gene belongs to the Rab family of monomeric GTPases, which are involved in the control of cellular membrane trafficking. The encoded protein plays a role in the targeted trafficking of integrins via its association with integrin alpha tails. As a consequence, the encoded protein is involved in the regulation of cell adhesion and migration. Expression of this gene is associated with a poor prognosis for glioma patients. This gene is downregulated by the tumor suppressor miR-200b, and miRNA-200b is itself downregulated in glioma tissues.
Applications:	WB, IHC, IF

Name of antibody:	RAB21
Immunogen:	Fusion protein of human RAB21
Full name:	RAB21, member RAS oncogene family (RAB21)
Synonyms:	KIAA0118
SwissProt:	Q9UL25
IHC positive control:	human prostate tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	24 kDa
WB Positive control:	HepG2, Hela, SVT2, A549, COS7, Jurkat, MDCK, PC12, MCF-7 cell lysates
WB Recommended dilution:	500-2000