

小鼠抗 METTL3 单克隆抗体

中文名称：小鼠抗 METTL3 单克隆抗体

英文名称：Anti-METTL3 mouse monoclonal antibody

别名：methyltransferase like 3; M6A; IME4; Spo8; MT-A70; hMETTL3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：METTL3

反应种属：Human, Rat, Mouse

标记物：Unconjugate

克隆类型：mouse monoclonal

技术规格

Background:	This gene encodes the 70 kDa subunit of MT-A which is part of N6-adenosine-methyltransferase . This enzyme is involved in the posttranscriptional methylation of internal adenosine residues in eukaryotic mRNAs, forming N6-methyladenosine.
Applications:	WB
Name of antibody:	METTL3
Immunogen:	Fusion protein of human METTL3
Full name:	methyltransferase like 3
Synonyms:	M6A; IME4; Spo8; MT-A70; hMETTL3

SwissProt:	Q86U44
WB Predicted band size:	64 KD
WB Positive control:	Raji, Hela, F9, RAW264.7, PC12, NIH/3T3 cell, Mouse kidney tissue lysates
WB Recommended dilution:	500-2000