

小鼠抗 MINOS1 单克隆抗体

中文名称: 小鼠抗 MINOS1 单克隆抗体

英文名称: Anti-MINOS1 mouse monoclonal antibody

别名: MIO10; Mic10; C1orf151

储存: 冷冻 (-20℃) 避光

抗原: MINOS1

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	Component of the MICOS complex, a large protein complex of the mitochondrial inner membrane that plays crucial roles in the maintenance of crista junctions, inner membrane architecture, and formation of contact sites to the outer membrane.
Applications:	WB, IHC
Name of antibody:	MINOS1
Immunogen:	Fusion protein of human MINOS1
Full name:	mitochondrial inner membrane organizing system 1 (MINOS1), transcript variant 1
Synonyms:	MIO10; Mic10; C1orf151

SwissProt:	Q5TGZ0
IHC positive control:	carcinoma of human lung tissue and adenocarcinoma of human endometrium tissue; human liver tissue and carcinoma of human liver tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	9 kDa
WB Positive control:	COS7 cell lysate
WB Recommended dilution:	500-2000