

小鼠抗 MLANA 单克隆抗体

中文名称: 小鼠抗 MLANA 单克隆抗体

英文名称: Anti-MLANA mouse monoclonal antibody

别名: MART-1; MART1

抗原: MLANA

储存: 冷冻 (-20°C) 避光

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	Involved in melanosome biogenesis by ensuring the stability of GPR143. Plays a vital role in the expression, stability, trafficking, and processing of melanocyte protein PMEL, which is critical to the formation of stage II melanosomes.
Applications:	WB, IHC, IF
Name of antibody:	MLANA
Immunogen:	Fusion protein of human MLANA
Full name:	melan-A (MLANA)
Synonyms:	MART-1; MART1
SwissProt:	Q16655

IHC positive control:	human pancreas tissue and human kidney tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	13 kDa
WB Positive control:	MDCK cell lysate
WB Recommended dilution:	500-2000