

## 小鼠抗 RDH14 单克隆抗体

中文名称: 小鼠抗 RDH14 单克隆抗体

英文名称: Anti-RDH14 mouse monoclonal antibody

别名: PAN2; SDR7C4

相关类别: 一抗

储存: 冷冻 (-20℃) 避光

宿主: Mouse

抗原: RDH14

反应种属: Human, Rat

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	Exhibits an oxidoreductive catalytic activity towards retinoids. Most efficient as an NADPH-dependent retinal reductase. Displays high activity toward 9-cis and all-trans-retinol. No steroid dehydrogenase activity detected.
<b>Applications:</b>	WB, FC
<b>Name of antibody:</b>	RDH14
<b>Immunogen:</b>	Fusion protein of human RDH14
<b>Full name:</b>	retinol dehydrogenase 14 (all-trans/9-cis/11-cis) (RDH14)
<b>Synonyms:</b>	PAN2; SDR7C4

<b>SwissProt:</b>	Q9HBH5
<b>WB Predicted band size:</b>	37 kDa
<b>WB Positive control:</b>	PC12 cell lysate
<b>WB Recommended dilution:</b>	200-1000