

小鼠抗 RELA 单克隆抗体

中文名称： 小鼠抗 RELA 单克隆抗体

英文名称： Anti-RELA mouse monoclonal antibody

别名： v-rel avian reticuloendotheliosis viral oncogene homolog A; p65; NFKB3

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： RELA

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	NF-kappa-B is a ubiquitous transcription factor involved in several biological processes. It is held in the cytoplasm in an inactive state by specific inhibitors. Upon degradation of the inhibitor, NF-kappa-B moves to the nucleus and activates transcription of specific genes. NF-kappa-B is composed of NFKB1 or NFKB2 bound to either REL, RELA, or RELB. The most abundant form of NF-kappa-B is NFKB1 complexed with the product of this gene, RELA. Four transcript variants encoding different isoforms have been found for this gene.
Applications:	WB, IHC

Name of antibody:	RELA
Immunogen:	Fusion protein of human RELA
Full name:	v-rel avian reticuloendotheliosis viral oncogene homolog A
Synonyms:	p65; NFKB3
SwissProt:	Q04206
IHC positive control:	Human tonsil and Human esophagus cancer; Human thyroid cancer and Human liver cancer
IHC Recommend dilution :	100-500
WB Predicted band size:	60 kDa
WB Positive control:	RAW264.7, 231, Hela cell lysates
WB Recommended dilution:	500-2000





