

小鼠抗 MS4A1 单克隆抗体

中文名称: 小鼠抗 MS4A1 单克隆抗体

英文名称: Anti-MS4A1 mouse monoclonal antibody

别名: B1; Bp35; CD20; CVID5; LEU-16; MS4A2; S7

抗原: MS4A1

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human , Dog , Monkey , Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:

This gene encodes a member of the membrane-spanning 4A gene family. Members of this nascent protein family are characterized by common structural features and similar intron/exon splice boundaries and display unique expression patterns among hematopoietic cells and nonlymphoid tissues. This gene encodes a B-lymphocyte surface molecule which plays a role in the development and differentiation of B-cells into plasma cells. This family member is localized to 11q12, among a cluster of family members. Alternative splicing of this gene results in two transcript variants which encode the same protein. [provi

	ded by RefSeq, Jul
Applications:	WB, IHC
Name of antibody:	MS4A1
Immunogen:	Fusion protein of human MS4A1
Full name:	membrane-spanning 4-domains, subfamily A, member 1 (MS4A1), transcript variant 3
Synonyms:	B1; Bp35; CD20; CVID5; LEU-16; MS4A2; S7
SwissProt:	P11836
IHC positive control:	human tonsil
IHC Recommend dilution:	30-150
WB Predicted band size:	33 kDa
WB Positive control:	SVT2, A549, COS7, MDCK cell lysates
WB Recommended dilution:	1000-5000