

小鼠抗 RHOJ 单克隆抗体

中文名称: 小鼠抗 RHOJ 单克隆抗体

英文名称: Anti-RHOJ mouse monoclonal antibody

别名: ARHJ; RASL7B; TC10B; TCL

抗原: RHOJ

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	This gene encodes one of the many small GTP-binding proteins in the Rho family shown to be associated with focal adhesions in endothelial cells (PMID: 21148427, 22103495). The encoded protein is activated by vascular endothelial growth factor and may regulate angiogenesis. [provided by RefSeq, Dec 2011].
Applications:	WB, IHC, IF
Name of antibody:	RHOJ
Immunogen:	Fusion protein of human RHOJ
Full name:	ras homolog family member J (RHOJ)
Synonyms:	ARHJ; RASL7B; TC10B; TCL

SwissProt:	Q9H4E5
IHC positive control:	human pancreas tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	24 kDa
WB Positive control:	Human breast tissue lysate
WB Recommended dilution:	500-2000