

小鼠抗 MSI1 单克隆抗体

中文名称： 小鼠抗 MSI1 单克隆抗体

英文名称： Anti-MSI1 mouse monoclonal antibody

别名： musashi RNA binding protein 1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： MSI1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:	This gene encodes a protein containing two conserved tandem RNA recognition motifs. Similar proteins in other species function as RNA-binding proteins and play central roles in posttranscriptional gene regulation. Expression of this gene has been correlated with the grade of the malignancy and proliferative activity in gliomas and melanomas. A pseudogene for this gene is located on chromosome 11q13.
Applications:	WB
Name of antibody:	MSI1
Immunogen:	Fusion protein of human MSI1

Full name:	musashi RNA binding protein 1
SwissProt:	O43347
WB Predicted band size:	39 kDa
WB Positive control:	U87MG, C6, IMR32 and PC12 cell, Mouse brain tissue lysates
WB Recommended dilution:	1000-5000

