

小鼠抗 RLIM 单克隆抗体

中文名称： 小鼠抗 RLIM 单克隆抗体

英文名称： Anti-RLIM mouse monoclonal antibody

别名： ring finger protein, LIM domain interacting; MRX61; RNF12; NY-REN-43

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： RLIM

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a RING-H2 zinc finger protein. It has been shown to be an E3 ubiquitin protein ligase that targets LIM domain binding 1 (LDB1/CLIM), and causes proteasome-dependent degradation of LDB1. This protein and LDB1 are co-repressors of LHX1/LIM-1, a homeodomain transcription factor. Multiple alternatively spliced variants, encoding the same protein, have been identified.
Applications:	WB, IHC
Name of antibody:	RLIM
Immunogen:	Fusion protein of human RLIM
Full name:	ring finger protein, LIM domain interacting

Synonyms:	MRX61; RNF12; NY-REN-43
SwissProt:	Q9NVW2
IHC positive control:	Human liver tissue
IHC Recommend dilution:	200-500
WB Predicted band size:	69 KD
WB Positive control:	HepG2, Hela, A549 cell lysates
WB Recommended dilution:	500-2000