

## 小鼠抗 MSMB 单克隆抗体

中文名称： 小鼠抗 MSMB 单克隆抗体

英文名称： Anti-MSMB mouse monoclonal antibody

别名： HPC13; IGBF; MSP; MSPB; PN44; PRPS; PSP; PSP-94; PSP57; PSP94

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Mouse

抗原： MSMB

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

#### Background:

The protein encoded by this gene is a member of the immunoglobulin binding factor family. It is synthesized by the epithelial cells of the prostate gland and secreted into the seminal plasma. This protein has inhibin-like activity. It may have a role as an autocrine paracrine factor in uterine, breast and other female reproductive tissues. The expression of the encoded protein is found to be decreased in prostate cancer. Two alternatively spliced transcript variants encoding different isoforms are described for this gene. The use of alternative polyadenylation sites has been found for this gene. [provided by RefSeq, Jul 2008].

<b>Applications:</b>	WB
<b>Name of antibody:</b>	MSMB
<b>Immunogen:</b>	Purified prostate specific antigen (MSMB) from human seminal fluid
<b>Full name:</b>	microseminoprotein, beta- (MSMB), transcript variant PSP94
<b>Synonyms:</b>	HPC13; IGBF; MSP; MSPB; PN44; PRPS; PSP; PSP-94; PSP57; PSP94
<b>SwissProt:</b>	P08118
<b>WB Predicted band size:</b>	13 kDa
<b>WB Positive control:</b>	MCF7 and Jurkat cell lysates
<b>WB Recommended dilution:</b>	200-1000