

## 兔抗 TOP2A(Ab-1343) 多克隆抗体

中文名称：兔抗 TOP2A(Ab-1343) 多克隆抗体

英文名称： Anti-TOP2A(Ab-1343) rabbit polyclonal antibody

别名： TOP2; TP2A

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： TOP2A(Ab-1343)

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	Control of topological states of DNA by transient breakage and subsequent rejoining of DNA strands. Topoisomerase II makes double-strand breaks. Essential during mitosis and meiosis for proper segregation of daughter chromosomes.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	TOP2A(Ab-1343)
<b>Immunogen:</b>	Synthesized non-phosphopeptide derived from human TOP2A around the phosphorylation site of threonine 1343 (E-K-T(p)-D-D).
<b>Full name:</b>	topoisomerase (DNA) II alpha 170kDa

<b>Synonyms :</b>	TOP2; TP2A
<b>SwissProt:</b>	P11388
<b>WB Predicted band size:</b>	174 kDa
<b>WB Positive control:</b>	HepG2 cells lysate
<b>WB Recommended dilution:</b>	500-3000

