

## 小鼠抗 SCRN1 单克隆抗体

中文名称: 小鼠抗 SCRN1 单克隆抗体

英文名称: Anti-SCRN1 mouse monoclonal antibody

别名: SES1

抗原: SCRN1

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human , Rat , Monkey

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	This gene likely encodes a member of the secernin family of proteins. A similar protein in rat functions in regulation of exocytosis in mast cells. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2009].
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	SCRN1
<b>Immunogen:</b>	Fusion protein of human SCRN1
<b>Full name:</b>	secernin 1 (SCRN1), transcript variant 2
<b>Synonyms:</b>	SES1
<b>SwissProt:</b>	Q12765

<b>IHC positive control:</b>	adenocarcinoma of human endometrium tissue and human lymph node tissue; human lymphoma tissue and adenocarcinoma of human colon tissue
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	46 kDa
<b>WB Positive control:</b>	HepG2, Hela, A549, COS7, Jurkat, PC12, MCF-7 cell lysates
<b>WB Recommended dilution:</b>	500-2000