

## 兔抗 TOP2B 多克隆抗体

中文名称：兔抗 TOP2B 多克隆抗体

英文名称： Anti-TOP2B rabbit polyclonal antibody

别名： TOPIIB; top2beta

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TOP2B

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a DNA topoisomerase, an enzyme that controls and alters the topologic states of DNA during transcription. This nuclear enzyme is involved in processes such as chromosome condensation, chromatid separation, and the relief of torsional stress that occurs during DNA transcription and replication. It catalyzes the transient breaking and rejoining of two strands of duplex DNA which allows the strands to pass through one another, thus altering the topology of DNA.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	TOP2B

<b>Immunogen:</b>	Synthetic peptide of human TOP2B
<b>Full name:</b>	topoisomerase (DNA) II beta 180kDa
<b>Synonyms :</b>	TOPIIB; top2beta
<b>SwissProt:</b>	Q02880
<b>ELISA Recommended dilution:</b>	2000-5000
<b>WB Predicted band size:</b>	183 kDa
<b>WB Positive control:</b>	Raw264.7, Jurkat, K562 and hepg2 cells
<b>WB Recommended dilution:</b>	500-2000

