

小鼠抗 SENP2 单克隆抗体

中文名称: 小鼠抗 SENP2 单克隆抗体

英文名称: Anti-SENP2 mouse monoclonal antibody

别名: AXAM2; SMT3IP2

抗原: SENP2

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human, Monkey

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	SUMO1 (UBL1; MIM 601912) is a small ubiquitin-like protein that can be covalently conjugated to other proteins. SENP2 is one of a group of enzymes that process newly synthesized SUMO1 into the conjugatable form and catalyze the deconjugation of SUMO1-containing species.[supplied by OMIM]. COMPLETENESS: complete on the 3' end.
Applications:	WB, IHC, IF, FC
Name of antibody:	SENP2
Immunogen:	Fusion protein of human SENP2
Full name:	SUMO1/sentrin/SMT3 specific peptidase 2 (SENP2)
Synonyms:	AXAM2; SMT3IP2

SwissProt:	Q9HC62
IHC positive control:	carcinoma of human lung tissue and adenocarcinoma of human ovary tissue; adenocarcinoma of human endometrium tissue and adenocarcinoma of human breast tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	68 kDa
WB Positive control:	HepG2, Hela, A549, COS7, Jurkat, MCF-7 cell lysates
WB Recommended dilution:	500-2000