

小鼠抗 SEPT1 单克隆抗体

中文名称: 小鼠抗 SEPT1 单克隆抗体

英文名称: Anti-SEPT1 mouse monoclonal antibody

别名: DIFF6; LARP; PNUTL3; SEP1

抗原: SEPT1

储存: 冷冻 (-20°C) 避光

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	This gene is a member of the septin family of GTPases. Members of this family are required for cytokinesis and the maintenance of cellular morphology. This gene encodes a protein that can form homo- and heterooligomeric filaments, and may contribute to the formation of neurofibrillary tangles in Alzheimer's disease. Alternatively spliced transcript variants have been found but the full-length nature of these variants has not been determined. [provided by RefSeq, Dec 2012].
Applications:	WB
Name of antibody:	SEPT1
Immunogen:	Fusion protein of human SEPT1

Full name:	septin 1 (SEPT1)
Synonyms:	DIFF6; LARP; PNUTL3; SEP1
SwissProt:	Q8WYJ6
WB Predicted band size:	42 kDa
WB Positive control:	Jurkat cell lysate
WB Recommended dilution:	200-1000