

## 小鼠抗 SIRT5 单克隆抗体

中文名称: 小鼠抗 SIRT5 单克隆抗体

英文名称: Anti-SIRT5 mouse monoclonal antibody

别 名: SIR2L5

相关类别: 一抗

储 存: 冷冻(-20℃) 避光

宿 主: Mouse

抗 原: SIRT5

反应种属: Human, Rat

标 记 物: Unconjugate

克隆类型: mouse monoclonal

## 技术规格

**Background:** 

This gene encodes a member of the sirtuin family of proteins, ho mologs to the yeast Sir2 protein. Members of the sirtuin family ar e characterized by a sirtuin core domain and grouped into four cl asses. The functions of human sirtuins have not yet been determined; however, yeast sirtuin proteins are known to regulate epigenetic gene silencing and suppress recombination of rDNA. Studies suggest that the human sirtuins may function as intracellular regulatory proteins with mono-ADP-ribosyltransferase activity. The protein encoded by this gene is included in class III of the sirtuin family. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq].



Applications:	WB, IHC, IF, FC
Name of antibody:	SIRT5
Immunogen:	Fusion protein of human SIRT5
Full name:	sirtuin 5 (SIRT5), transcript variant 1
Synonyms:	SIR2L5
SwissProt:	Q9NXA8
IHC positive control:	human kidney tissue
IHC Recommend dilut ion:	30-150
WB Predicted band si ze:	34 kDa
WB Positive control:	HepG2, Hela, COS7, PC12, MCF-7 cell lysates
WB Recommended dil ution:	200-1000