

小鼠抗 SIRT6 单克隆抗体

中文名称： 小鼠抗 SIRT6 单克隆抗体

英文名称： Anti-SIRT6 mouse monoclonal antibody

别名： sirtuin 6; SIR2L6

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： SIRT6

反应种属： Human, Mouse, Rat, Monkey


标记物： Unconjugate



克隆类型： Mouse Monoclonal

技术规格

Background:	This gene encodes a member of the sirtuin family of NAD-dependent enzymes that are implicated in cellular stress resistance, genomic stability, aging and energy homeostasis. The encoded protein is localized to the nucleus, exhibits ADP-ribosyl transferase and histone deacetylase activities, and plays a role in DNA repair, maintenance of telomeric chromatin, inflammation, lipid and glucose metabolism. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	WB, IF, IP
Name of antibody:	SIRT6
Immunogen:	Fusion protein of human SIRT6

Full name:	sirtuin 6
Synonyms:	SIR2L6
SwissProt:	Q8N6T7
WB Predicted band size:	39 kDa
WB Positive control:	Jurkat, MCF7, HCT116, K562, COS7, C2C12, 3T3 and HeLa cell lysates
WB Recommended dilution:	500-2000
IF positive control:	HeLa cells





mlbio 梅联生物
Good elisakit producers