

## 小鼠抗 SMAP2 单克隆抗体

中文名称： 小鼠抗 SMAP2 单克隆抗体

英文名称： Anti-SMAP2 mouse monoclonal antibody

别名： small ArfGAP2; SMAP1L

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： SMAP2

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	GTPase activating protein that acts on ARF1. Can also activate ARF6 (in vitro). May play a role in clathrin-dependent retrograde transport from early endosomes to the trans-Golgi network (By similarity).
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	SMAP2
<b>Immunogen:</b>	Fusion protein of human SMAP2
<b>Full name:</b>	small ArfGAP2
<b>Synonyms:</b>	SMAP1L
<b>SwissProt:</b>	Q8WU79
<b>IHC positive control:</b>	Human lymphoma tissue

<b>IHC Recommend dilution:</b>	200-500
<b>WB Predicted band size:</b>	47 KD
<b>WB Positive control:</b>	A431 and U251 cell lysates
<b>WB Recommended dilution:</b>	500-2000