

小鼠抗 SMARCA4 单克隆抗体

- 中文名称: 小鼠抗 SMARCA4 单克隆抗体
- 英文名称: Anti-SMARCA4 mouse monoclonal antibody


别 名: SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 4; BRG1; CSS4; SNF2; SWI2; MRD16; RT

- 相关类别: 一抗
- 储 存: 冷冻(-20℃)
- 宿 主: Mouse
- 抗 原: SMARCA4
- 反应种属: Human, Mouse
- 标记物: Unconjugate
- 克隆类型: Mouse Monoclonal


技术规格



	bdoid tumor predisposition syndrome type 2. Multiple transcript variants encoding different isoforms have been found for this g ene.
Applications:	WB, IF, IP
Name of antibody:	SMARCA4
Immunogen:	Fusion protein of human SMARCA4
Full name:	SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 4
Synonyms:	BRG1; CSS4; SNF2; SWI2; MRD16; RTPS2; BAF190; SNF2L4; SNF2 LB; hSNF2b; BAF190A
SwissProt:	P51532
WB Predicted band size :	185 kDa
WB Positive control:	K562, F9 and Hela NE cell lysates
WB Recommended dilu tion:	1000-5000
IF positive control:	Hela cells
IF Recommend dilution:	25-50







SMARCA4

