

## 兔抗 TRAF7 多克隆抗体

中文名称：兔抗 TRAF7 多克隆抗体

英文名称： Anti-TRAF7 rabbit polyclonal antibody

别名： RFWD1, RNF119

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TRAF7

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Tumor necrosis factor (TNF; see MIM 191160) receptor-associated factors, such as TRAF7, are signal transducers for members of the TNF receptor superfamily (see MIM 191190). TRAFs are composed of an N-terminal cysteine/histidine-rich region containing zinc RING and/or zinc finger motifs; a coiled-coil (leucine zipper) motif; and a homologous region that defines the TRAF family, the TRAF domain, which is involved in self-association and receptor binding.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	TRAF7
<b>Immunogen:</b>	Synthetic peptide of human TRAF7

<b>Full name:</b>	TNF receptor-associated factor 7, E3 ubiquitin protein ligase
<b>Synonyms :</b>	RFWD1, RNF119
<b>SwissProt:</b>	Q6Q0C0
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human colon cancer and Human gastric cancer
<b>IHC Recommend dilution:</b>	15-50
<b>WB Predicted band size:</b>	67 kDa
<b>WB Positive control:</b>	Jurkat cells
<b>WB Recommended dilution:</b>	200-1000



