

小鼠抗 SRPRB 单克隆抗体

中文名称: 小鼠抗 SRPRB 单克隆抗体

英文名称: Anti-SRPRB mouse monoclonal antibody

别名: APMCF1

储存: 冷冻 (-20℃) 避光

抗原: SRPRB

宿主: Mouse

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	The protein encoded by this gene has similarity to mouse protein in which is a subunit of the signal recognition particle receptor (SR). This subunit is a transmembrane GTPase belonging to the GTPase superfamily. It anchors alpha subunit, a peripheral membrane GTPase, to the ER membrane. SR is required for the cotranslational targeting of both secretory and membrane proteins to the ER membrane. [provided by RefSeq, Jul 2008].
Applications:	WB, IHC
Name of antibody:	SRPRB
Immunogen:	Fusion protein of human SRPRB
Full name:	SRP receptor beta subunit (SRPRB)

Synonyms:	APMCF1
SwissProt:	Q9Y5M8
IHC positive control:	human lymphoma tissue and human tonsil tissue
IHC Recommend dilution:	100-500
WB Predicted band size:	30 kDa
WB Positive control:	Hela, Jurkat, MCF-7 cell lysates
WB Recommended dilution:	500-2000