

小鼠抗 SSB 单克隆抗体

中文名称： 小鼠抗 SSB 单克隆抗体

英文名称： Anti-SSB mouse monoclonal antibody

别名： Sjogren syndrome antigen B; La; LARP3; La/SSB

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： SSB

反应种属： Human

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:	The protein encoded by this gene is involved in diverse aspects of RNA metabolism, including binding and protecting poly(U) termini of nascent RNA polymerase III transcripts from exonuclease digestion, processing 5' and 3' ends of pre-tRNA precursors, acting as an RNA chaperone, and binding viral RNAs associated with hepatitis C virus. Autoantibodies reacting with this protein are found in the sera of patients with Sjogren syndrome and systemic lupus erythematosus. Alternative promoter usage results in two different transcript variants which encode the same protein.
Applications:	WB

Name of antibody:	SSB
Immunogen:	Fusion protein of human SSB
Full name:	Sjogren syndrome antigen B
Synonyms:	La; LARP3; La/SSB
SwissProt:	P05455
WB Predicted band size :	47 kDa
WB Positive control:	MCF7, A549, Jurkat and Hela cell lysates
WB Recommended dilution:	1000-5000

