

小鼠抗 STAT6 单克隆抗体

中文名称: 小鼠抗 STAT6 单克隆抗体

英文名称: Anti-STAT6 mouse monoclonal antibody

别 名: signal transducer and activator of transcription 6; STAT6B; STAT6C; D12S1644;

IL-4-STAT

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Mouse

抗 原: STAT6

反应种属: Human

标 记 物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:

The protein encoded by this gene is a member of the STAT f amily of transcription factors. In response to cytokines and gr owth factors, STAT family members are phosphorylated by the receptor associated kinases, and then form homo- or heterodi mers that translocate to the cell nucleus where they act as tr anscription activators. This protein plays a central role in exert ing IL4 mediated biological responses. It is found to induce t he expression of BCL2L1/BCL-X(L), which is responsible for the anti-apoptotic activity of IL4. Knockout studies in mice sugges ted the roles of this gene in differentiation of T helper 2 (Th2



) cells, expression of cell surface markers, and class switch of immunoglobulins. Alternative splicing results in multiple transc ript variants.
Applications:	ELISA, WB
Name of antibody:	STAT6
Immunogen:	Fusion protein of human STAT6
Full name:	signal transducer and activator of transcription 6
Synonyms:	STAT6B; STAT6C; D12S1644; IL-4-STAT
SwissProt:	P42226
WB Predicted band size:	94 kDa
WB Positive control:	HeLa whole cell lysate
WB Recommended dilutio n:	1000-5000

