

小鼠抗 STAU1 单克隆抗体

- 中文名称: 小鼠抗 STAU1 单克隆抗体
- 英文名称: Anti-STAU1 mouse monoclonal antibody
- 别名: staufen double-stranded RNA binding protein 1; STAU; PPP1R150
- 相关类别: 一抗
- 储存: 冷冻(-20℃)
- 宿 主: Mouse
- 抗 原: STAU1
- 反应种属: Human
- 标记物: Unconjugate
- 克隆类型: mouse monoclonal

技术规格

	Staufen is a member of the family of double-stranded RNA (ds RNA)-binding proteins involved in the transport and/or localizat ion of mRNAs to different subcellular compartments and/or org anelles. These proteins are characterized by the presence of mu Itiple dsRNA-binding domains which are required to bind RNAs
Background:	having double-stranded secondary structures. The human homo logue of staufen encoded by STAU, in addition contains a micr otubule- binding domain similar to that of microtubule-associat ed protein 1B, and binds tubulin. The STAU gene product has been shown to be present in the cytoplasm in association with the rough endoplasmic reticulum (RER), implicating this protein



	in the transport of mRNA via the microtubule network to the R ER, the site of translation. Five transcript variants resulting from alternative splicing of STAU gene and encoding three isoforms have been described. Three of these variants encode the same isoform, however, differ in their 5'UTR.
Applications:	WB
Name of antibody:	STAU1
Immunogen:	Fusion protein of human STAU1
Full name:	staufen double-stranded RNA binding protein 1
Synonyms:	STAU; PPP1R150
SwissProt:	O95793
WB Predicted band size:	63 KD
WB Positive control:	Jurkat, K562, HepG2, Hela cell lysates
WB Recommended diluti on:	500-2000