

## 小鼠抗 STUB1 单克隆抗体

中文名称： 小鼠抗 STUB1 单克隆抗体

英文名称： Anti-STUB1 mouse monoclonal antibody

别名： STIP1 homology and U-box containing protein 1; CHIP; UBOX1; SCAR16; HSPABP2; NY-CO-7; SDCCAG7

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： STUB1

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein containing tetratricopeptide repeat and a U-box that functions as a ubiquitin ligase/cochaperone. The encoded protein binds to and ubiquitinates shock cognate 71 kDa protein (Hspa8) and DNA polymerase beta (Polb), among other targets. Mutations in this gene cause spinocerebellar ataxia, autosomal recessive 16. Alternative splicing results in multiple transcript variants. There is a pseudogene for this gene on chromosome 2.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	STUB1

<b>Immunogen:</b>	Fusion protein of human STUB1
<b>Full name:</b>	STIP1 homology and U-box containing protein 1
<b>Synonyms:</b>	CHIP; UBOX1; SCAR16; HSPABP2; NY-CO-7; SDCCAG7
<b>SwissProt:</b>	Q9UNE7
<b>WB Predicted band size:</b>	35 kDa
<b>WB Positive control:</b>	HeLa whole cell lysate
<b>WB Recommended dilution:</b>	500-2000

