

## 小鼠抗 STX3 单克隆抗体

中文名称:	小鼠抗 STX3 单克隆抗体
英文名称:	Anti-STX3 mouse monoclonal antibody
别名:	syntaxin 3; STX3A
相关类别:	一抗
储存:	冷冻 (-20℃)
宿主:	Mouse
抗原:	STX3
反应种属:	Human, Rat, Mouse
标记物:	Unconjugate
克隆类型:	mouse monoclonal

### 技术规格

<b>Background:</b>	The gene is a member of the syntaxin family. The encoded protein is targeted to the apical membrane of epithelial cells where it forms clusters and is important in establishing and maintaining polarity necessary for protein trafficking involving vesicle fusion and exocytosis. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	STX3
<b>Immunogen:</b>	Fusion protein of human STX3
<b>Full name:</b>	syntaxin 3
<b>Synonyms:</b>	STX3A

<b>SwissProt:</b>	Q13277
<b>WB Predicted band size:</b>	33 KD
<b>WB Positive control:</b>	HepG2, 293T, RAW264.7, PC12, NIH/3T3 cell, Human kidney tissue lysates
<b>WB Recommended dilution:</b>	500-2000