

## 小鼠抗 SULT1A3 单克隆抗体

中文名称: 小鼠抗 SULT1A3 单克隆抗体

英文名称: Anti-SULT1A3 mouse monoclonal antibody

别 名: sulfotransferase family 1A member 3; STM; HAST; HAST3; M-PST; ST1A3; ST1A5; TL-PST; ST1A3/ST1A4

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Mouse

抗原: SULT1A3

反应种属: Human

标记物: Unconjugate

克隆类型: mouse monoclonal

## 技术规格

**Background:** 

Sulfotransferase enzymes catalyze the sulfate conjugation of many hormones, neurotransmitters, drugs, and xenobiotic co mpounds. These cytosolic enzymes are different in their tissu e distributions and substrate specificities. The gene structure (number and length of exons) is similar among family memb ers. This gene encodes a phenol sulfotransferase with thermo labile enzyme activity. Four sulfotransferase genes are located on the p arm of chromosome 16; this gene and SULT1A4 ar ose from a segmental duplication. This gene is the most cen tromeric of the four sulfotransferase genes. Read-through tra



	nscription exists between this gene and the upstream SLX1A (SLX1 structure-specific endonuclease subunit homolog A) ge ne that encodes a protein containing GIY-YIG domains.
Applications:	WB
Name of antibody:	SULT1A3
Immunogen:	Fusion protein of human SULT1A3
Full name:	sulfotransferase family 1A member 3
Synonyms:	STM; HAST; HAST3; M-PST; ST1A3; ST1A5; TL-PST; ST1A3/ST1 A4
SwissProt:	P0DMN0
WB Predicted band size:	34 KD
WB Positive control:	HepG2 cell lysate
WB Recommended dilutio n:	500-2000