

小鼠抗 SWAP70 单克隆抗体

中文名称： 小鼠抗 SWAP70 单克隆抗体

英文名称： Anti-SWAP70 mouse monoclonal antibody

别名： SWAP switching B-cell complex subunit 70; HSPC321; SWAP-70

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： SWAP70

反应种属： Human, Mouse, Rat, Monkey

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:	Phosphatidylinositol 3,4,5-trisphosphate-dependent guanine nucleotide exchange factor (GEF) which, independently of RAS, transduces signals from tyrosine kinases to RAC. It also mediates signaling of membrane ruffling. Regulates the actin cytoskeleton as an effector or adapter protein in response to agonist stimulated phosphatidylinositol (3,4)-bisphosphate production and cell protrusion (By similarity).
Applications:	WB, IP
Name of antibody:	SWAP70
Immunogen:	Fusion protein of human SWAP70

Full name:	SWAP switching B-cell complex subunit 70
Synonyms:	HSPC321; SWAP-70
SwissProt:	Q9UH65
WB Predicted band size:	69 kDa
WB Positive control:	COS7, 293T, C6, 3T3 and Hela cell lysates
WB Recommended dilution:	1000-5000



