

小鼠抗 SYP 单克隆抗体

中文名称: 小鼠抗 SYP 单克隆抗体

英文名称: Anti-SYP mouse monoclonal antibody

别名: MRX96; MRXSYP

抗原: SYP

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	This gene encodes an integral membrane protein of small synaptic vesicles in brain and endocrine cells. The protein also binds cholesterol and is thought to direct targeting of vesicle-associated membrane protein 2 (synaptobrevin) to intracellular compartments. Mutations in this gene are associated with X-linked mental retardation (XLMR). [provided by RefSeq, Aug
Applications:	WB, IF
Name of antibody:	SYP
Immunogen:	Fusion protein of human SYP
Full name:	synaptophysin (SYP)
Synonyms:	MRX96; MRXSYP

SwissProt:	P08247
WB Predicted band size:	34 kDa
WB Positive control:	Hela cell lysate
WB Recommended dilution:	200-1000