

## 小鼠抗 TBCC 单克隆抗体

中文名称:	小鼠抗 TBCC 单克隆抗体
英文名称:	Anti-TBCC mouse monoclonal antibody
别名:	CFC
储存:	冷冻 (-20℃) 避光
抗原:	TBCC
宿主:	Mouse
反应种属:	Human , Rat , Monkey , Mouse
相关类别:	一抗
标记物:	Unconjugate
克隆类型:	mouse monoclonal

### 技术规格

<b>Background:</b>	Cofactor C is one of four proteins (cofactors A, D, E, and C) involved in the pathway leading to correctly folded beta-tubulin from folding intermediates. Cofactors A and D are believed to play a role in capturing and stabilizing beta-tubulin intermediates in a quasi-native confirmation. Cofactor E binds to the cofactor D/beta-tubulin complex; interaction with cofactor C then causes the release of beta-tubulin polypeptides that are committed to the native state. [provided by RefSeq].
<b>Applications:</b>	WB, IHC, IF
<b>Name of antibody:</b>	TBCC

<b>Immunogen:</b>	Fusion protein of human TBCC
<b>Full name:</b>	tubulin folding cofactor C (TBCC)
<b>Synonyms:</b>	CFC
<b>SwissProt:</b>	Q15814
<b>IHC positive control:</b>	adenocarcinoma of human ovary tissue and human pancreas tissue; adenocarcinoma of human endometrium tissue and a denocarcinoma of human colon tissue
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	39 kDa
<b>WB Positive control:</b>	HepG2, Hela, SVT2, A549, COS7, PC12, MCF-7 cell lysates
<b>WB Recommended dilution:</b>	500-2000