

小鼠抗 TDO2 单克隆抗体

中文名称: 小鼠抗 TDO2 单克隆抗体

英文名称: Anti-TDO2 mouse monoclonal antibody

别名: TDO; TO; TPH2; TRPO

储存: 冷冻 (-20℃) 避光

抗原: TDO2

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	Tryptophan 2,3-dioxygenase (EC 1.13.11.11) plays a role in catalyzing the first and rat-limiting step in the kynurenine pathway, the major pathway of tryptophan metabolism.[supplied by OMIM].
Applications:	WB, IHC
Name of antibody:	TDO2
Immunogen:	Fusion protein of human TDO2
Full name:	tryptophan 2,3-dioxygenase (TDO2)
Synonyms:	TDO; TO; TPH2; TRPO
SwissProt:	P48775
IHC positive control:	adenocarcinoma of human breast tissue

IHC Recommend dilution:	30-150
WB Predicted band size:	48 kDa
WB Positive control:	Hela, A549 cell lysates
WB Recommended dilution:	500-2000