

## 小鼠抗 TEKT2 单克隆抗体

中文名称： 小鼠抗 TEKT2 单克隆抗体

英文名称： Anti-TEKT2 mouse monoclonal antibody

别名： tektin 2; TEKTB1; TEKTIN-T; h-tektin-t

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： TEKT2

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	This gene product belongs to the tektin family of proteins. Tektins comprise a family of filament-forming proteins that are coassembled with tubulins to form ciliary and flagellar microtubules. This gene is expressed in the testis and its protein is localized to the flagella of the sperms, indicating that it may play a role in spermatogenesis.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	TEKT2
<b>Immunogen:</b>	Fusion protein of human TEKT2
<b>Full name:</b>	tektin 2
<b>Synonyms:</b>	TEKTB1; TEKTIN-T; h-tektin-t



<b>SwissProt:</b>	Q9UIF3
<b>WB Predicted band size:</b>	50 KD
<b>WB Positive control:</b>	Hela and 293T cell lysates
<b>WB Recommended dilution:</b>	500-2000